

Chronos

Refraction System

Guided Binocular Refraction.

Streamline your workflow and delegate refraction with Chronos



**COMPACT, RELIABLE REFRACTION
SYSTEM** that combines binocular
autorefractometry and keratometry with
binocular subjective testing and
visual acuity.

Chronos

Refraction System



Streamline your workflow and delegate refraction with Chronos

Chronos is a multifunctional, space-saving instrument
that optimises your workflow.

O V E R V I E W



DELEGATE

SightPilot™ is a guided refraction system that simplifies exams and facilitates delegation.



STREAMLINE YOUR WORKFLOW

Divert straightforward patients through Chronos, reserving time for more complex cases. Customise the Chronos exam depending on whether it is pre-operative, post-operative, conventional refraction, etc.



SAVE SPACE

Chronos, as an all-in-one platform with a small footprint, avoids the need to factor in the patient-chart distance in the room layout, saving space and boosting cost efficiency, and providing flexibility on where refraction takes place.



SAVE TIME

Chronos saves time by optimising the workflow, eliminating the time needed to clean and move between devices.

DISCOVER HOW YOU CAN STREAMLINE YOUR WORKFLOW AND DELEGATE REFRACTION WITH CHRONOS



CHRONOS- ENHANCED WORKFLOW




Time Saving



Maximising
Resource Allocation

Accuracy of refraction is paramount whether you are refracting pre or post-surgery, or undertaking routine refraction. Chronos allows you to delegate refraction without compromise, for straightforward patients.



→ Current workflow

→ Additional pre-operative workflow option by adding Chronos

→ Additional post-operative workflow option by adding Chronos

1

TOPCON

Patient Details

Previous spectacles

Input prescription of previous spectacles

Right	Sphere + -3.25	Cylinder + -0.75	Axis + 92	ADD + 1.00
Left	Sphere + -3.50	Cylinder + -1.00	Axis + 98	ADD + 1.00

Save

Back

Quit

Reset

Enter the patient information and import/enter the previous spectacle prescription to begin the refraction¹.

Figure 1 displays four screenshots of the AR application interface. The top-left screenshot shows the 'Instructume' screen, which displays a 3D model of a microscope and a yellow line indicating the correct viewing angle. The top-right screenshot shows the 'Instructions' screen, which displays the same 3D model with a yellow arrow indicating the correct viewing angle. The bottom-left screenshot shows the 'Vortex' screen, which displays a 3D model of a vortex with a yellow line indicating the correct viewing angle. The bottom-right screenshot shows the 'Vortex' screen, which displays a 3D model of a vortex with a yellow line indicating the correct viewing angle.

SightPilot™ provides step-by-step instructions to position the patient and then automatically aligns the optics to complete the objective refraction.

Figure 1 displays four screenshots of the TAPCON software interface, illustrating the calibration process:

- Top Left:** The "New Gun Compass" screen shows a grid of colored dots (red, green, blue) used for initial alignment. It includes "R" and "L" buttons for rotation and a "Test" button.
- Top Right:** The "Gun Angle Adjustment" screen shows a target with a cluster of dots. It includes "R" and "L" buttons for rotation and a "Test" button.
- Bottom Left:** The "Gun Power Adjustment" screen shows a target with a cluster of dots. It includes "R" and "L" buttons for rotation and a "Test" button.
- Bottom Right:** The "Visual Audit" screen shows a target with a cluster of dots. It includes "R" and "L" buttons for rotation and a "Test" button.

SightPilot™ walks the operator through a variety of subjective refraction tests including visual acuity charts, red-green comparison, cylinder adjustment, binocular balancing and near addition charts. On-screen prompts enable quick input of patient response to advance to the next step in the process.

The screenshot shows the TOPCON website with a purple header. The main content area is titled 'Results' and displays data for a 'Significant reduction (subjective reduction)'.

Significant reduction (subjective reduction)

	2014-2015	2016-2017	2018	2019	1st	2nd
Rate	0.50	-1.50	1.75	2.50	15	15
Flow	0.75	-0.75	45	7.50	15	15
Flow distance	60.0	56.5	The distance of 1st floor distance (1st floor distance) is 100m.			

Additional test information

Subjective reduction

	2014-2015	2016	1st
Rate	0.50	-1.50	1.75
Flow	0.75	-0.75	45
Flow distance	60.0		

Stratified

	2014-2015	2016-2017	2018	2019	1st
Rate					25

Buttons at the bottom: Back, Manual control, Finish exam.

When the refraction is complete, the results are displayed on screen and may be printed or sent to the patient's EHR file.

1. Previous prescription can be entered manually or imported directly from your Topcon lensmeter (SOLOS [Manufactured by VISIA IMAGING S.R.L.] or CL-300)



Chronos


Refraction System

With Chronos **automated binocular refraction system**, spend more time on what matters most, your patients.




Chronos binocular refraction technology

Chronos measures both autorefractometry and subjective refraction under binocular viewing conditions, for a more natural, comfortable visual experience. Binocular refraction has been shown to provide better control of accommodation for objective and subjective end points.



Reduce alignment errors with Chronos auto-alignment

Chronos uses Topcon's 3D stereo camera technology to optimise alignment throughout testing, pioneered in Topcon's automated OCTs and retinal cameras.



Cutting-edge moving lens system

Chronos incorporates a cutting-edge moving lens system enabling rapid and smooth changes in spherical lens power. This provides a more comfortable visual experience for the patient.




Ease of test distance adjustment

A combination of the built-in θ movement mechanism of the head and the lens movement, adjusts the convergence angle and adjusts the different testing distances.

TEST DISTANCE

Far-/Near-point test distance can be set

25cm - 609.6cm



Far


Intermediate

Near


*Conventional refraction set-up, for comparison

Accuracy with Chronos


A combination of the super luminescent diode (SLD) ring, rotary prism technology and binocular objective refraction, provides stable measurements, including in patients with some media opacification.



Measured simultaneously



Patented Rotary Prism Technology



Conventional ring

Chronos ring

Specifications

Objective measurement

Refraction measurement range	Spherical refractive power	-25D ~ +22D ^{*1,2}
	Cylindrical refractive power	-10D ~ 0D ^{*1,2}
	Cylinder axis angle	1° ~ 180°
Corneal curvature measurement range	Corneal curvature radius	5.00mm ~ 10.00mm
	Corneal refractive power	67.50D ~ 33.75D (Conversion value when the corneal refractive ratio is 1.3375)
	Corneal principal meridian angle	1° ~ 180°
Minimum measurement unit	Spherical/cylindrical refractive power	0.12D
	Cylinder axis angle	1°
	Corneal curvature radius	0.01mm
	Corneal refractive power	0.12D
	Corneal principal meridian angle	1°
Display of measured value	Displayed on the control screen of the operation controller	
Minimum measurable pupil diameter	Φ2.0mm	
PD measurement range	50mm ~ 80mm	
Minimum PD measurement unit	0.5mm	

Subjective measurement

Refraction measurement range	Spherical refractive power	-18.00D ≤ Equivalent spherical power ≤ +18.00D ^{*3}
	Cylindrical refractive power	-8.00D ≤ Cylindrical refractive power (Cylindrical power) ≤ 0.00D ^{*4}
	<i>All conditions stated on the right must be met^{*5}</i>	
	Cylinder axis angle	1° ~ 180°
	Horizontal prism (one eye movable range)	±15.0Δ ^{*6}
Minimum measurement unit	Vertical prism (one eye movable range)	±2.5Δ
	Spherical/ADD refractive power	0.25D
	Cylindrical refractive power	0.25D
	Cylinder axis angle	1°
Test distance	Prism refractive power	0.1Δ
	Far-/Near-point test distance can be set between 25cm and 6.096m	
Visual acuity measurement range ^{*7}	0.05 ~ 1.6 (decimal notation)	
Chart	Visual acuity test chart; spherical power correction test chart, astigmatism test chart and binocular function test chart	
Background luminance	155±15cd/m ²	
Display of measured value	Displayed on the screen of the operation controller	
Record of measured value	Printing by thermal printer/external printer, data output	
Measuring head movement	Right-and-left direction	Inside 9mm to Outside 12.5mm
	Up-and-down direction	Down 15mm to Up 15mm
	Back-and-forth direction	Forward: 20mm - Backward: 20mm
Measuring head rotary angle	Convergence 17.5° to Divergence 8.5° (Eyeball torsion axis centre)	

Other Specifications

Dimensions and Weight	Main unit	Dimensions: 510-540mm (H) × 671-766mm (W) × 278-357mm (D) Weight: 31.2 kg
	Power supply unit	Dimensions: 276mm (H) x 117mm (W) x 197mm (D) Weight: 3.5 kg
Electric Rating	Source voltage	AC100 ~ 240V
	Frequency	50 ~ 60Hz
	Power input	160VA

*1 The dioptric powers are indicated with reference wavelength $\lambda_0 = 587.56 \text{ nm}$

*2 Spherical refractive power + Cylindrical refractive power ≤ +22D or Spherical refractive power + Cylindrical refractive power ≥ -25D

*3 The conversion value with "VD=12mm" is described here.

*4 The conversion value with "VD=3mm" is described here.

*5 The value described here is the maximum value. The measurement range is smaller according to the test distance setting for executing a test or the setting conditions of VD during measurement.

*6 The value described here is the maximum value. The measurable range is smaller according to the combination of the patient's PD and the test distance.

*7 0.1 ~ 1.6 complies with ISO 10938. ETDRS chart using Landolt Ring (visual acuity 0.25 ~ 1.6) complies with ANSI Z80.21.

Optional Accessories



SNI-1

SightPilot NAVi™ Controller SNI-1 for REFRACTION SYSTEM Chronos

SightPilot NAVi™ Patient-Guided Refraction with Voice Prompts and a Simple Controller SNI-1.

* The external speaker or headphones are not included and must be purchased locally. Please check with your local distributor for details.



CRX-1000

Chronos Chinrest

For more stable measurement, an optional chinrest attachment is available.

IMPORTANT In order to obtain the best results with this instrument, please be sure to review all user instructions prior to operation. Not all products, services, or offers are available in all markets. Contact your local distributor for country-specific information and availability.

CLASS 1 LASER PRODUCT
(IEC60825-1:2014)
PRODUIT LASER DE CLASSE 1
(CEI60825-1:2014)

TOPCON INSTRUMENTS (MALAYSIA) SDN. BHD.
(Regional Office for Topcon Healthcare Southeast Asia)
Unit 2, 4, Jalan Pensyarah U1/28, Hicom-glenmarie Industrial Park,
40150 Shah Alam, Selangor, MALAYSIA
Phone: +603-766 16260 Fax: +603-766 16261
Email: mys_tim_marketing_sm@topcon.com
www.topconhealthcare.my

TOPCON SINGAPORE MEDICAL PTE. LTD.
100G Pasir Panjang Road, #02-18, Interlocal Centre,
SINGAPORE 118523
Phone: +65-68720606 Fax: +65-67736150
E-mail: med.sales.sg@topcon.com
www.topconhealthcare.sg

TOPCON INSTRUMENTS (THAILAND) CO., LTD.
77/162 Sinnsathorn Tower, 37th Floor, Krungthoburi Rd.,
Klongtongsa, Klongsarn, Bangkok 10600, THAILAND
Phone: +66-02-440-1152 Fax: +66-02-440-1158
Email: tha_medical@topcon.com
www.eyecare.topcon.co.th

MEHRA EYETECH PRIVATE LIMITED
801 B Wing, Lotus Corporate Park, Graham Firth Steel Compound
Goregaon (East) Mumbai 400063 Maharashtra, INDIA
Phone: +91-22-61285455
E-mail: sales@mehraeyetech.in
www.topconhealthcare.in

TOPCON (BEIJING) MEDICAL TECHNOLOGY CO., LTD.
Room 2808, Tower C, JinChangAn Building, No.82, Middle Section of
East 4th Ring Road, Chaoyang District, Beijing 100124, P.R. CHINA
Phone: +86-10-8794-5176
E-mail: cn_marketing@topcon.com
www.topcon-china.net/

TOPCON KOREA MEDICAL CO., LTD.
2F YK Building, 205, Dogok-ro, Gangnam-gu,
Seoul, Republic of Korea
Phone: +82-2-6959-7947
E-mail: tkm@topcon.com
www.topconhealthcare.kr/



TOPCON CORPORATION

75-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8580, JAPAN.
Phone: +81-(0)3-3558-2522/2502
Fax: +81-(0)3-3965-6898
www.topconhealthcare.jp

TOPCON Healthcare